


Best Seller

2~8°C Pharmacy Refrigerator

Intelligent Control & Constant Temperature
 Multiple Alarms Bring More Security



YC-55L



Temp. Uniformity



Pulling Down to 5°C



Temp. Recovery After Door Opened for 1 Min

**Typical freezer data based on internal testing with freezer set point at 5°C and ambient temperature at 22°C.*

Features



Temperature Control System
 Display temp. in 0.1°C increments



Compact and Fit
 Fits under standard 34" high ADA compliant countertop



Hydrocarbon Refrigerant
 Environment-friendly Natural Gas



CE Certificate





Double-pane Glass Door
 Automatic Electric Heating Function



Alarms System

Structure





2~8°C Pharmacy Refrigerator

| | |
|-------------------------------|-------------------------------|
| Model | YC-55L |
| Capacity | 55L/1.94Cu.Ft |
| Internal Size (W*D*H) mm/Inch | 444*440*404/17.48*17.32*15.90 |
| External Size (W*D*H) mm/Inch | 540*560*632/21.20*22.00*24.90 |
| Package Size (W*D*H) mm/Inch | 575*617*682/22.64*24.30*26.85 |
| NW/GW (kgs) | 35/41 |

Performance

| | |
|---------------------|-----------------|
| Temperature Range | 2~8°C |
| Ambient Temperature | 16~32°C |
| Cooling Performance | 5°C |
| Climate Class | N |
| Controller | Microprocessor |
| Display | Digital Display |

Refrigeration

| | |
|----------------------|--------------------|
| Compressor | 1pc |
| Cooling Method | Forcedair Cooling |
| Defrost Mode | Automatical |
| Refrigerant | R600a |
| Insulation Thickness | D:55, L/R:48, B:50 |

Construction

| | |
|---------------------|---|
| External Material | PCM |
| Inner Material | Sprayed Aluminum Plate |
| Shelves | 2(Plastic Coated Steel Wire Shelves) |
| Access Test Port | 1pc. Ø23mm |
| Casters | 2+(2 Leveling Feet) |
| Data Logging/Export | Record every 10 mins/Max.12 months via USB port |
| Interface | Remote Alarm Contact, RS485 |
| Backup Battery | Yes |

Alarms

| | |
|-------------|--|
| Temperature | High/Low Temperature, High Ambient Temperature |
| Electrical | Power Failure, Low Battery, Sensor Error |
| System | Data Communication Failure, Door Ajar, Data Logger Failure |

Electrical

| | | | |
|--------------------|-----------|--------------|--------------|
| Power Supply(V/Hz) | 110V/60Hz | 220~240/50Hz | 220~240/60Hz |
| Rated Current(A) | 1.16A | 0.9A | 0.56A |

*Zhongke Meiling reserves the right to change products and specifications without prior notice.